

ISLAND SCIENTIFIC LTD

Reply to:

Unit 9
Ventnor Industrial Estate
Station Road
Ventnor
Isle of Wight
PO38 1DX
England



*Vacuum Equipment, Used and Reconditioned for Technology,
Industry, Research and Education*

Tel No: (0) 1983 855822
Fax No: (0) 1983 852146

E-Mail Address: sales@island-scientific.co.uk
Web Site: <http://www.island-scientific.co.uk>



Make: Leybold

Model: Eco Dry L

Pumping speed: 48CMH

Ultimate Vacuum: 3 x 10⁻² m/bar

Inlet flange: NW40

Outlet flange: NW25

Cooling: Air

Voltage: 415V

Motor Power: 1500W

Overall Dimensions: 755 mm long x 572 mm high x 220 mm wide

Weight: 95 Kg

Condition: Used but in Good Condition

Accessories Available: