

# ISLAND SCIENTIFIC LTD



*Vacuum Equipment, Used and Reconditioned for Technology,  
Industry, Research and Education*

*Reply to:*

Unit 9  
Ventnor Industrial Estate  
Station Road  
Ventnor  
Isle of Wight  
PO38 1DX  
England

Tel No: (0) 1983 855822

Fax No: (0) 1983 852146

E-Mail Address: [sales@island-scientific.co.uk](mailto:sales@island-scientific.co.uk)

Web Site: <http://www.island-scientific.co.uk>



***Make: Welch***

***Model: 2580C-50 Dry Pump***

*Displacement: 4.9 m/hr*

*Vacuum: 12 m/bar*

*Pressure: N/A*

*Inlet connection: 10 mm Tubing*

*Outlet connection: Hose tail*

*Motor power: 0.25Kw, 240V, 1 Phase*

*Type of oil: N/A*

*Weight: 7.5Kg*

*Length: 282 mm*

*Width: 165 mm*

*Height: 231mm*

*Accessories available:*