

ISLAND SCIENTIFIC LTD



*Vacuum Equipment, Used and Reconditioned for Technology,
Industry, Research and Education*

Reply to:

Unit 9
Ventnor Industrial Estate
Station Road
Ventnor
Isle of Wight
PO38 1DX
England

Tel No: (0) 1983 855822

Fax No: (0) 1983 852146

E-Mail Address: sales@island-scientific.co.uk

Web Site: <http://www.island-scientific.co.uk>



Make: Edwards

Model: PRCT10-C

Part No. D021-53-000

Type: Pirani

Range: 5 to 10⁻² m/bar

Connection flange: Cone Tube

Height: 155 mm

Width: 40 mm

Depth: 26 mm

Weight: 0.2 kg

Operating Temperature Range: 5-60°C

Suitable gauge controller: Pirani 12 or 161 model supplied separately