

ISLAND SCIENTIFIC LTD

Reply to:

Unit 4
Ventnor Industrial Estate
Station Road
Ventnor
Isle of Wight
PO38 1DX
England

Tel No: (0) 1983 855822


Fax No: (0) 1983 852146

E-Mail Address: sales@island-scientific.co.uk

Web Site: <http://www.island-scientific.co.uk>



*Vacuum Equipment, Used and Reconditioned for Technology,
Industry, Research and Education*



Make: Pfeiffer

Model: WKD500

Displacement: 490CMH

Inlet connection: ISO 100

Outlet connection: ISO 100

Type of oil: H/C

Minimum size backing pump: 40 CMH

Motor power: 1.5 kW, 380-415V, 3 Phase

Cooling: Air

Weight: 115 Kg

Length: 790 mm

Width: 324 mm

Height: 285 mm