

ISLAND SCIENTIFIC LTD



*Vacuum Equipment, Used and Reconditioned for Technology,
Industry, Research and Education*

Reply to:

Unit 9
Ventnor Industrial Estate
Station Road
Ventnor
Isle of Wight
PO38 1DX
England

Tel No: (0) 1983 855822

Fax No: (0) 1983 852146

E-Mail Address: sales@island-scientific.co.uk

Web Site: <http://www.island-scientific.co.uk>



Make: Edwards

Model: XDS10C

Pumping speed: 11.3 CMH

Ultimate Vacuum: 6.8×10^{-1} m/bar

Inlet flange: KF25

Outlet flange: KF25

Cooling: Air

Voltage: 240V

Motor Power: .3 Kw

Overall Dimensions: 315 mm long x 288 mm high x 249 mm wide

Weight: 24.5 Kg

Condition: Reconditioned

Accessories Available: