

ISLAND SCIENTIFIC LTD

Reply to:

Unit 9
Ventnor Industrial Estate
Station Road
Ventnor
Isle of Wight
PO38 1DX
England

Tel No: (0) 1983 855822
Fax No: (0) 1983 852146

E-Mail Address: sales@island-scientific.co.uk

Web Site: <http://www.island-scientific.co.uk>



*Vacuum Equipment, Used and Reconditioned for Technology,
Industry, Research and Education*



Make: Edwards

Model: EXC300

Part No. D396-14-000

Turbo Pump Controller for EXT Turbo Pumps

Electrical Supply Voltage: 240v

Maximum Input Power (Without Bake out Band): 550w

Maximum Output Power: 290w

Ambient Temperature Range: 0-40°C

Weight: 4.6

Width: 213mm

Height: 132mm

Length: 296mm