

# ISLAND SCIENTIFIC LTD

Reply to:

Unit 9  
Ventnor Industrial Estate  
Station Road  
Ventnor  
Isle of Wight  
PO38 1DX  
England

Tel No: (0) 1983 855822

Fax No: (0) 1983 852146

E-Mail Address: [sales@island-scientific.co.uk](mailto:sales@island-scientific.co.uk)

Web Site: <http://www.island-scientific.co.uk>



*Vacuum Equipment, Used and Reconditioned for  
Technology, Industry, Research and Education*



Un-serviced Unit

***Make: Varian***

***Model: HS20***

*Displacement: 17,500 lit/sec*

*Inlet flange: ISO630*

*Backing connection: ISO100*

*Cooling: Water*

*Fluid capacity: 4.7 litres*

*Type of fluid: Variable*

*Recommended minimum backing pump: 170CMH*

*Heater power: 12kw 415V 3Phase*

*Diameter:*

*Height: 1219mm*

*Weight: Kg*